

Stability Studies to Support PET Drug Applications

Danny Bingham
Director PET cGMP Compliance
Triad Isotopes

Objective

My objective is to provide information to help you decide the best stability study practices for your submission. We will discuss the following:

- Review the purpose for stability studies
- Review what to include in stability studies
- Process development/control relationship to stability studies
- Vial size and orientation considerations
- Expiration time
- Elevated temperatures and other considerations

What are Stability Studies?

Stability studies are used to confirm the drug product is:

- stable under the defined storage conditions
- stable at the highest radioactive concentration
- stable at expiration time

The entire batch volume is stored in the intended container/closure for the duration of the study.

Stability Studies

Stability studies address concerns due to:

- Radiation-related radiolysis
- Chemical changes
 - Interaction between drug product components and vial/stopper
 - Ability of container/closure to protect against contamination (container/closure integrity)

Analytical methods must distinguish degradation products and impurities

- Process development and process controls – impurity profile
 - Analytical method development/validation

Stability Studies – New Submission

Three stability batches (qualification batches) are provided in the initial submission

- Should be performed on three consecutive business days – adequate planning
- Why not all in one day?
 - ALARA concerns
 - Logistics – difficult to perform 3 stability batches in one day (especially for F-18 products)
 - You are not involving parts of the quality system that may impact the chemical profile. For example:
 - Daily reagent/diluent preparation
 - Daily chemistry module set-up/cleaning
 - Formulation dilution
 - Daily mobile phase preparation
 - Cleaning

Stability Studies – New Submission

Stability program is described in your submission

- Includes a summary paragraph
- Table summarizing QC tests and stability indicating tests for the product
 - Summary of the stability batch data.
- Post-approval commitment
 - Single stability study performed annually

How do I develop the program?

- What are your goals?
 - Derived from process development
- Once these are defined then the protocol can be written.

What do I need to include in stability studies?

The stability protocol typically consists of the following:

1. Summary/purpose, Schedule, Procedure
2. Full QC at T_0 (as defined in your submission)
3. Perform stability indicating tests at defined time points:
 - Radiochemical identity and purity
 - Radionuclidic purity
 - Appearance
 - pH
 - Stabilizer/preservative effectiveness
 - Chemical purity
 - Specific activity (if necessary)
4. Summary report
 - Includes complete and legible batch records

Best Practices

Best Practices

- Appearance
- pH
- Radiochemical purity and identity
- Stabilizer/preservative effectiveness
- Chemical purity
- Specific activity (if necessary)
- Radionuclidic purity
- Endotoxin testing
- Sterility testing


Best practices for stability study QC testing begins with process development.

- Establishes which QC and stability indicating tests to perform (or not perform)


If you claim it then prove it.

- Most CMC review questions answered with process development data...

Best Practices



Best Practices – Precursor Stability



Best Practices – Process Development

Take the time to ensure process development is adequate:

- There is value in “justifying the ruggedness” of a process with a certain number of batches before finalizing your CMC section.
 - Especially for complex syntheses
 - Finalize process controls/stability program
 - The stability batches must be performed using the process submitted in the application
 - Do not want to be in the position of amending the CMC section of an application during review
 - This will require the stability study/qualification batches to be repeated
- Provides support for expiration time and time points.

Best Practices – Expiration

What about expiration time?

- Decision is based on stability goals - process development
 - A good starting point for new products is 6 half-lives

Which time points should I use for stability testing?

- Well studied products (for example FDG, NaF, ammonia) have used two points: EOS (t_0) and t_{expiry}
 - Two points have been used successfully for IND drug products
- An additional time point may be needed
 - For example: t_0 , t_6 , t_{12}
 - Decision is based on process development/process controls

What About Vial Orientation?

Should vials be stored upside down or right side up?

Stability studies should be performed with inverted configuration

- Formulation is in contact with all surfaces of the container/closure
 - Interaction between drug product components and vial/stopper
 - Container/closure integrity
 - Specified in the PET drug application guidance “Sample Formats – Chemistry, Manufacturing, and Controls (CMC) Section” for FDG, NaF, ammonia.

Upright configuration may be acceptable.

- Justify in the application

What About Vial Size?

Do you need to perform stability studies on each vial size used in PET drug manufacturing?

The PET Q&A Guidance for PET Drugs (Dec 2012) addresses this (Q68):

- Recommends that the largest vial size be used
- If you manufacture in 30 mL and 50 mL vials, then use 50 mL vials
 - May need to test the 30 mL vial if the headspace oxygen-to-surface ratio differs significantly
 - If you manufacture into 30 mL vials and want to use 50 mL vials, then new stability studies will be required (in triplicate)
- Will require a prior-approval supplement for a change in the container closure.

Other Considerations

Elevated temperature study – necessary?

- Likely not needed
 - PET Q&A Guidance for PET Drugs (Dec 2012) Q69: not required for commonly used PET drugs (e.g. FDG, NaF, ammonia)
 - New drug products: may need to consider elevated temperatures.
 - Especially if shipped long distances
 - Use process development data to justify your position.
- Other stability considerations
 - Product is photosensitive

Discuss with FDA review division – For example, end of Phase 2

Other Considerations

Stability studies to be performed in triplicate for:

- Change in strength
- Change in stabilizer/preservative content
- Change in final product container/closure
- Change in storage conditions
- Change in expiration time
- Addition of a precursor manufacturer
- Addition of a target material manufacturer (for example, [^{18}O]water)

Submitted as an annual report, CBE or prior approval supplement as appropriate.

Documentation

Lack of good documentation practices is a frequent cause of stability study issues during internal reviews

Some errors have resulted in the rejection of a stability batch.

- A thorough review at T_0 and T_{expiry} will save time correcting mistakes
 - or starting over
- Second person review critical

Risk during FDA review...

Documentation

Consider the following good documentation practices:

1. Clear – Must be legible
2. Traceable – (who, what, when, where, why)
3. Accurate – Avoid transcription errors
4. Complete – The record must be complete with all attachments
5. Reviewed – Timely and detailed review of the batch record is critical
 - Errors detected can be immediately corrected
 - Procrastination can lead to exasperation
 - Annotated properly – Single line through/date/initial and explanation

Summary

Best practices are derived from process development and can be used to support/justify these stability study goals:

- Impurity profile - analytical methods must distinguish degradation products/impurities
- Vial orientation
- Expiration time
- Elevated temperatures/photosensitivity/other

Stability program - defines t_0 and t_{expiry} stability testing to perform, or not to perform

- Consistent with goals
- Stability protocol derived from the stability program
- If you claim it, prove it

Good documentation practices are important

- Don't jeopardize your study (or submission) with errors
- A 30 minute detailed review can save hours of time

Thank You!